

# DAT Facts™

ASTM D6228

<b>Analyte:</b>	<b>CAS #</b>	<b>ppb</b>
Hydrogen sulfide	773-06-4	5
Carbonyl sulfide	463-58-1	5
Methyl mercaptan	74-93-6	5
Ethyl mercaptan	75-08-1	5
Dimethyl sulfide	75-17-3	5
Carbon disulfide	75-15-0	5
2-Propyl mercaptan	75-33-2	5
Thiirane	420-12-2	5
tert-Butyl mercaptan	75-66-1	5
n-Propyl mercaptan	107-03-9	5
Ethyl methyl sulfide	624-89-5	5
Thiophene	110-02-1	5
iso-Butyl mercaptan	513-44-0	5
Diethyl sulfide	110-81-6	5
n-Butyl mercaptan	109-79-5	5
Dimethyl disulfide	75-18-3	5
3-Methyl thiophene	616-44-4	5
Tetrahydrothiophene	110-01-0	5
2-Ethyl thiophene	872-55-9	5
2,5-Dimethyl thiophene	638-02-8	5
Diethyl disulfide	110-81-6	5
Allyl sulfide	592-88-1	5